



IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!
Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.
Should you give this device to another person, it is vital that you also pass on these instructions for use.

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GB - ENGLISH

- Using the Lexie:**
- Warm/cold pillow for the body.
 - Warm: relaxes and reduces pain in the abdominal area during menstruation.
 - Neck and shoulder wrap.
 - Pain relief for lower back pain
 - Relaxation: the warmth and weight will help you relax.
 - During pregnancy: as a heat pad for the relief of ligament pain and lower back pain (LBP)
 - Also great for soothing aching muscles.
 - Ideal as a cold compress for muscle aches, swelling, insect bites,

- sunburn, headaches, fever, summer heat and hot flushes.
- As a heat pillow: provides soothing heat to relieve tension, stress and pain.
- It is suitable for a wide range of uses and conforms well to the neck, shoulders, back, head, feet and abdomen.
- The Lexie is fragrance-free.

- WARNING:**
Always follow these instructions and keep them in a safe place. Always read the instructions before use.
- Always use with the cover
 - The pillow should not be used by anyone who is hypersensitive or insensitive to heat and other very vulnerable people who are unable to respond to overheating.
 - Do not expose the pillow to excessive heat from other heat sources. The product is intended for indoor use.
 - ~~✂~~ Make sure it does not get wet! The pillow can no longer be used if it gets wet.
 - To ensure a long service life, it is very important to set the microwave correctly when it is heated.
 - ~~🚫~~ Not recommended for children under 4 years of age or for use on individuals with medical conditions, who are disabled, asleep or insensitive to heat (those unable to respond to overheating).
 - ~~👁~~ Always use under supervision.
 - Do not use on people who are asleep
 - The pillow should only be used on a covered part of the body. This is NOT a toy. Parental supervision is recommended.
 - Do not heat the pillow in a regular oven or a convection oven so those methods of heating are more difficult to control.
 - Do not let the pillow come into contact with metal while it is being heated. For example, touch the grill, the wall or the bottom of the microwave. This will damage the pillow.
 - Always check the temperature of the pillow with the inside of your wrist before use.
 - Avoid contact with sensitive skin.
 - Do not use if there are any signs of skin irritation.
 - Do not use if there are any symptoms of inflammation.
 - The pillow itself is NOT washable. However, the cover can be washed up to 40 degrees.
 - Do not use the pillow if it is worn or has overheated (the seeds will char and give off an unpleasant baking smell, or the pillow cover will show scorch marks).
 - The cotton fabric will be affected and irreparably damaged if the padding gets too hot and the seeds retain too much heat. This can usually be seen by the brown spots and small holes that appear in the fabric. The pillow is then unusable. Moreover, the burnt smell is also unpleasant.
 - Caution: Never leave the pillow unattended while heating it in an oven or microwave. There is a risk of fire if the heating process is not carried out properly.
 - KipKep cannot be held responsible for any damage caused by using the Lexie improperly or by failing to follow the instructions for using the pillow.
 - Any serious incident that occurs when using the pillow should be reported to the manufacturer (KipKep B.V.) and the competent authority of the country where the user/patient resides.
 - Always follow your microwave manufacturer's instructions to be sure of the power level for each setting on your microwave.
 - Composition. Filling: linseed | Pillow: 100% cotton | Cover: varies, see label in cover. Always use the pillow with the original KipKep cover.
 - Check the pillow regularly for signs of any wear or damage. If there are any signs that the pillow has been used incorrectly or is not working, contact the supplier before using it again.

- User Instructions KipKep Lexie Heating Pillow**
Can be used as a heat pad and as a cold compress. The shape of the Lexie makes it ideal for use on the stomach, lower back and neck. The weight is distributed across the pillow. The pillow should feel lukewarm, never hot.
- ~~✂~~ Always remove the grill from your microwave!
 - Do not exceed the recommended heating time and do not reheat the pillow until it has cooled down to room temperature.
 - Always heat the pillow on the rotating glass turntable. Preferably on a clean microwave-proof (ceramic or glass) plate. Always consult the user manual of your microwave.
 - Make sure your microwave is free of any leftover food. Leftover food becomes very hot when heated and may cause scorch marks in the fabric of the pillow.
 - Let the pillow cool for 1 minute before placing it on somebody's body. After a minute, always check the temperature on the inside of your wrist or in the crook of your arm (max 45 C).
 - Always heat the pillow without a cover.
 - Always use the pillow in its cover. Placed the heated or cooled pillow in the cover.

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- Never use the pillow directly on bare skin.
 - Wash the cover before first use. To clean it, you can wash the outer cover at up to 40 C degrees.
 - The pillow stays comfortably warm for up to half an hour and can be reheated. Allow the pillow to completely cool down between uses.
- Heating:**
- (Fig.1) Remove the pillow from the cover.
 - Always remove the grill from the microwave before heating the pillow and make sure the pillow is on a clean surface.
 - Never place it on a metal surface and make sure it cannot touch the microwave walls during heating. Carefully place the pillow in the microwave and ensure that nothing obstructs the movement of the plate.
 - Fold the pillow in half; spread the seeds evenly. If your microwave does not have a glass turntable, place the pillow on a microwave-safe plate/bowl and rotate it halfway through the heating time to make sure that the product heats up evenly. Shake the pillow back and forth once in between.
 - Make sure that features such as grilling or browning (if available on the microwave) are switched off.
 - (Fig. 1) Heat the pillow without the cover in a microwave at 600 Watt for 90 seconds. If necessary, increase the time by 15 seconds until the desired temperature has been reached. Never heat the pillow for more than 2 minutes!
 - Do not leave the pillow unattended while it is being heated.
 - Carefully remove the pillow from the microwave to avoid burning yourself.
 - A microwave works by slowly spreading heat from the inside to the outside. This means it will take a while for the pillow to feel warm after being removed from the microwave.
 - OVEN: heat for 15 minutes at 100 degrees in a clean baking dish. The pillow should feel lukewarm, not too hot to touch. Overheating may cause fire or scorching and may damage the pillow and the surroundings.
 - (Fig. 2) Wait 1 minute after heating it and check the temperature of the pillow whether it has cooled sufficiently before use.
 - Always check the temperature as described above.
 - (Fig. 2) Put the pillow back in the cover.
 - Do not leave the pillow in the microwave for too long to avoid overheating!
 - Always use under supervision.

- Cooling:** When used as a cold compress.
- Cool the pillow with cover in a clean freezer bag in the fridge or freezer.
 - Place the pillow in the freezer for at least 1 hour.
 - When refrigerated, always place the pillow in a clean freezer bag to avoid contact with food and/or ice.
- Storage:**
- When the Lexie is not in use, store it in a dry and cool place.
 - If you store the pillow, let it cool down completely before storing it.
 - Do not store in direct sunlight.

- CE-Markering:**
This product complies with the requirements and other relevant provisions of the current European CE Directives, Medical Devices, Class 1

- Warranty:**
The local and applicable legislation governing the warranty applies to your purchase. KipKep B.V. cannot accept any warranty claims in the following cases: damage or defects caused by use; operation or handling of the product in ways that are not consistent with normal personal or household use; nor (attempted) repairs by a party other than KipKep B.V.
- IMPORTANT:** Changes or modifications that have not been explicitly approved by the party responsible for handling these may lead to the loss of the right of the user to use this product.

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